

Biodesign®

Dit biologiske valg

Biodesign®
ADVANCED TISSUE REPAIR

cookbiotech.eu

| | |
|--|-----------|
| Introduktion | 3 |
| Vores historie | 3 |
| Et fundament af løbende forbedringer | 4 |
| Teknologi | 5 |
| SIS-teknologi | 5 |
| Kompleks sammensætning | 6 |
| Biotech Process™ | 7 |
| Sikker | 12 |
| Vævsspecifik omdannelse | 12 |
| Pålidelig | 14 |
| Undersøgt og påvist | 14 |
| Vigtig klinisk evidens | 15 |
| Produkter | 16 |
| Biodesign® 4-Layer Tissue Graft | 16 |
| Biodesign® Dural Graft | 16 |
| Biodesign® Duraplasty Graft | 16 |
| Biodesign® Fistula Plug Set | 17 |
| Biodesign® Hernia Graft | 17 |
| Biodesign® Hiatal Hernia Graft | 17 |
| Biodesign® Otologic Repair Graft | 18 |
| Biodesign® Rectopexy Graft | 18 |
| Litteratur | 19 |

Vores historie

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Vigtig klinisk evidens

Abdominal wall reconstruction with components separation and mesh reinforcement in complex hernia repair

Nockolds CL, Hodde JP, Rooney PS. Abdominal wall reconstruction with components separation and mesh reinforcement in complex hernia repair. *BMC Surg*. 2014;14:25.

The use of porcine small intestinal submucosa as a prosthetic material for laparoscopic hernia repair in infected and potentially contaminated fields: Long-term follow-up

Franklin ME Jr, Treviño JM, Portillo G, Vela I, Glass JL, González JJ. The use of porcine small intestinal submucosa as a prosthetic material for laparoscopic hernia repair in infected and potentially contaminated fields: Long-term follow-up. *Surg Endosc*. 2008;22(9):1941-1946.

Biologic prosthesis to prevent recurrence after laparoscopic paraesophageal hernia repair: long-term follow-up from a multicenter, prospective, randomized trial

Oelschläger BK, Pellegrini CA, Hunter JG, et al. Biologic prosthesis to prevent recurrence after laparoscopic paraesophageal hernia repair: long-term follow-up from a multicenter, prospective, randomized trial. *J Am Coll Surg*. 2011;213(4):461-468.

Porcine small intestinal submucosa (SIS) myringoplasty in children: a randomized controlled study


D'Eredita RD. Porcine small intestinal submucosa (SIS) myringoplasty in children: a randomized controlled study. *Int J Pediatr Otorhinolaryngol*. 2015;79(7):1085-1089.

Ventral herniorrhaphy: experience with two different biosynthetic mesh materials, Surgisis and Alloderm

Gupta A, Zahriya K, Mullens PL, Salmassi S, Keshishian A. Ventral herniorrhaphy: experience with two different biosynthetic mesh materials, Surgisis and Alloderm. *Hernia*. 2006;10(5):419-425.

Porcine small intestine submucosal graft for endoscopic skull base reconstruction

Illing E, Chaaban MR, Riley KO, Woodworth BA. Porcine small intestine submucosal graft for endoscopic skull base reconstruction. *Int Forum Allergy Rhinol*. 2013;3(11):928-932.

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





Illustrationer af Lisa Clark

1. Internal Cook Biotech document: D00278582.
2. Clause KC, Barker TH. Extracellular matrix signaling in morphogenesis and repair. *Curr Opin Biotechnol.* 2013;24(5):830-833.
3. Daley WP, Peters SB, Larsen M. Extracellular matrix dynamics in development and regenerative medicine. *J Cell Sci.* 2008;121(Pt 3):255-264.
4. Hodde J, Janis A, Ernst D, Zopf D, Sherman D, Johnson C. Effects of sterilization on an extracellular matrix scaffold: Part I. Composition and matrix architecture. *J Mater Sci Mater Med.* 2007;18(4):537-543.
5. Nihsen ES, Johnson CE, Hiles MC. Bioactivity of small intestinal submucosa and oxidized regenerated cellulose/collagen. *Adv Skin Wound Care.* 2008;21(10):479-486.
6. Hodde J, Janis A, Hiles M. Effects of sterilization on an extracellular matrix scaffold: Part II. Bioactivity and matrix interaction. *J Mater Sci Mater Med.* 2007;18(4):545-550.
7. Turner NJ, Badylak SF. Biologic scaffolds for musculotendinous tissue repair. *Eur Cell Mater.* 2013;25:130-143.
8. Franklin ME Jr, Trevino JM, Portillo G, Vela I, Glass JL, Gonzalez JJ. The use of porcine small intestinal submucosa as a prosthetic material for laparoscopic hernia repair in infected and potentially contaminated field: Long-term follow-up. *Surg Endosc.* 2008;22(9):1941-1946.
9. Stelly M, Stelly TC. Histology of CorMatrix bioscaffold 5 years after pericardial closure. *Ann Thorac Surg.* 2013;96(5):e127-e129.
10. Badylak S, Kokini K, Tullius B, Whitson B. Strength over time of a resorbable bioscaffold for body wall repair in a dog model. *J Surg Res.* 2001;99(2):282-287.
11. Badylak SF, Park K, Peppas N, McCabe G, Yoder M. Marrow-derived cells populate scaffolds composed of xenogeneic extracellular matrix. *Exp Hematol.* 2001;29(11):1310-1318.
12. Hodde J. Extracellular matrix as a bioactive material for soft tissue reconstruction. *ANZ J Surg.* 2006;76(12):1096-1100.
13. Badylak SF. The extracellular matrix as a scaffold for tissue reconstruction. *Semin Cell Dev Biol.* 2002;13(5):377-383.
14. Frantz C, Stewart KM, Weaver VM. The extracellular matrix at a glance. *J Cell Sci.* 2010;123(24):4195-4200.
15. Hiles M, Record Ritchie RD, Altizer AM. Are biologic grafts effective for hernia repair?: a systematic review of the literature. *Surg Innov.* 2009;16(1):26-37.
16. Gupta A, Zahriya K, Mullens PL, Salmassi S, Keshishian A. Ventral herniorrhaphy: experience with two different biosynthetic mesh materials, Surgisis and Alloderm. *Hernia.* 2006;10(5):419-425.
17. Kissane NA, Itani KMF. A decade of ventral incisional hernia repairs with biologic acellular dermal matrix: What have we learned? *Plast Reconstr Surg.* 2012;130(5 Suppl 2):194S-202S.
18. Novitsky YW, Rosen MJ. The biology of biologics: basic science and clinical concepts. *Plast Reconstr Surg.* 2012;130(5 Suppl 2):9S-17S.
19. Internal Cook Biotech document: 00-107.
20. Internal Cook Biotech document: 02-063.
21. Swinehart IT, Badylak SF. Extracellular matrix bioscaffolds in tissue remodeling and morphogenesis. *Dev Dyn.* 2016;245(3):351-360.
22. Raeder RH, Badylak SF, Sheehan C, Kallakury B, Metzger DW. Natural anti-galactose alpha1,3 galactose antibodies delay, but do not prevent the acceptance of extracellular matrix xenografts. *Transpl Immunol.* 2002;10(1):15-24.
23. Bach FH, Ferran C, Hechenleitner P, et al. Accommodation of vascularized xenografts: expression of "protective genes" by donor endothelial cells in a host Th2 cytokine environment. *Nat Med.* 1997;3(2):196-204.

24. Allman AJ, McPherson TB, Badylak SF, et al. Xenogeneic extracellular matrix grafts elicit a TH2-restricted immune response. *Transplantation*. 2001;71(11):1631-1640.
25. Allman AJ, McPherson TB, Merrill LC, Badylak SF, Metzger DW. The Th2-restricted immune response to xenogeneic small intestinal submucosa does not influence systemic protective immunity to viral and bacterial pathogens. *Tissue Eng*. 2002;8(1):53-62.
26. Badylak SF, Valentin JE, Ravindra AK, McCabe GP, Stewart-Akers AM. Macrophage phenotype as a determinant of biologic scaffold remodeling. *Tissue Eng Part A*. 2008;14(11):1835-1842.
27. Martinez FO, Sica A, Mantovani A, Locati M. Macrophage activation and polarization. *Front Biosci*. 2008;13:453-461.
28. Bochicchio GV, Jain A, McGonigal K, et al. Biologic vs. synthetic inguinal hernia repair: 1-year results of a randomized double-blinded trial. *J Am Coll Surg*. 2014;218(4):751-757.
29. Ansaloni L, Catena F, Coccolini F, Gazzotti F, D'Alessandro L, Pinna AD. Inguinal hernia repair with porcine small intestine submucosa: 3-year follow-up results of a randomized controlled trial of Lichtenstein's repair with polypropylene mesh vs. Surgisis Inguinal Hernia Matrix. *Am J Surg*. 2009;198(3):303-312.
30. Smart NJ, Pathak S, Boorman P, Daniels IR. Synthetic or biological mesh use in laparoscopic ventral mesh rectopexy - a systematic review. *Colorectal Dis*. 2013;15(6):650-654.
31. Albayati S, Morgan MJ, Turner CE. Laparoscopic ventral rectopexy for rectal prolapse and rectal intussusception using a biological mesh. *Colorectal Dis*. 2017;19(8):857-862.
32. Oelschlager BK, Pellegrini CA, Hunter JG, et al. Biologic prosthesis to prevent recurrence after laparoscopic paraesophageal hernia repair: long-term follow-up from a multicenter, prospective, randomized trial. *J Am Coll Surg*. 2011;213(4):461-468.

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