COOK
MEDICAL

CE 0123

EN
3

**Biodesign®
Tissue Graft**

Instructions for Use

CS
4

**Biodesign®
Tkáňový graft**

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DA
6

**Biodesign®
vænstransplantat**

Brugsanvisning

DE
7

**Biodesign®
Gewebeimplantat**

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EL
9

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
BIODESIGN® TISSUE GRAFT

DEVICE DESCRIPTION:

The Biodesign® 4-Layer Tissue Graft is a dried multi-layered small intestinal submucosa (SIS) sheet. The graft may be implanted in all patient populations to reinforce soft tissue and can be cut to size to accommodate the patient's anatomy. Over the course of several months, the patient's own healing mechanisms incorporate the SIS into the adjacent tissues until it can no longer be detected, resulting in new, remodeled tissue.^{1,3} Long-term follow-up cystoscopy (36-62 weeks post-operation) has shown that SIS-grafted areas appear as normal tissue, as confirmed by histological examination.³ The 4-Layer Tissue Graft is provided sterile for single use only.

INTENDED USE

The Biodesign 4-Layer Tissue Graft is intended for implantation to reinforce soft tissue. It is indicated for use in chordee correction, Peyronie's disease treatment, and urethral repair. The graft is supplied sterile in peel-open packages and is intended for one-time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

This symbol means the following: Magnetic Resonance Safe

4-LAYER TISSUE GRAFT This symbol means the following: 4-Layer Tissue Graft

This product is intended for use by medical professionals trained in urologic surgical procedures.

CONTRAINDICATIONS

The graft is derived from a porcine source and should not be used for patients with known sensitivity to porcine materials.

PRECAUTIONS

- The graft is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to graft failure and/or transmission of disease.
- **Do not resterilize.** Discard all open and unused portions.
- The graft is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard the graft if mishandling has caused possible damage or contamination, or if the graft is past its expiration date.
- Ensure that the graft is hydrated prior to suturing or stapling.
- Patients undergoing radiation therapy may not experience normal wound healing.
- **Graft performance has not been evaluated with suture spacing greater than 2 mm.**
- **Ensure that all layers of the graft are secured when suturing or stapling.**

POTENTIAL COMPLICATIONS

Complications that can occur with the use of any surgical graft material may include:

- Abscess
- Adhesion
- Allergic Reaction
- Bleeding
- Delayed or failed incorporation
- Erectile dysfunction
- Erosion
- Extrusion
- Fever
- Fistula Formation
- Hematoma
- Induration
- Infection
- Inflammation
- Migration
- Pain
- Persistence or recurrence of tissue defect
- Seroma formation
- Tissue trauma

If any of the above conditions occur and cannot be resolved, device removal may be considered.

STORAGE

Store the graft in a clean, dry location at room temperature.

STERILIZATION

The graft has been sterilized with ethylene oxide. Do not resterilize.

INSTRUCTIONS FOR USE

Required Materials

- Sterile basin
- Sterile smooth forceps
- Hydration fluid: room temperature sterile saline or sterile lactated Ringer's solution

NOTE: Handle the graft using aseptic technique, minimizing contact with latex gloves.

1. Remove the graft inner pouch from its outer pouch, and place the inner pouch onto the sterile field.
2. Open the inner pouch and aseptically remove the graft with the sterile smooth forceps.
3. Place the graft into the sterile basin on the sterile field. (Multiple grafts may be hydrated simultaneously in the same basin.)
4. Add the hydration fluid to the basin.
5. Hydrate the graft in the fluid until the desired handling characteristics are achieved. A hydration time of greater than 1 minute is not required. Do not hydrate the graft for longer than 5 minutes.
6. Prepare the implant site using standard surgical techniques.
7. Trim the graft to fit the site of the implant, providing a small allowance for overlap.

NOTE: The graft may be trimmed prior to hydration to fit the patient's anatomy. Be sure to hydrate the graft prior to suturing it into place. See step 5.

8. After hydration, suture or staple the graft into place, avoiding excess tension.

NOTE: Surgical experience indicates that suturing or stapling graft sheets with close tissue approximation results in better outcomes.

9. Complete the surgical procedure.
10. Discard any unused portions of the graft according to institutional guidelines for medical waste.

REFERENCES

1. Mantovani F, Tondelli E, Cozzi G, et al. Reconstructive urethroplasty using porcine acellular matrix (SIS): evolution of the grafting technique and results of 10-year experience. *Urologia*. 2011;78(2):92-97.
2. Palminteri E, Berdondini E, Fusco F, De Nunzio C, Salonia A. Long-term results of small intestinal submucosa graft in bulbar urethral reconstruction. *Urology*. 2012;79(3):695-701.
3. Xu YM, Fu Q, Sa YL, Zhang J, Song LJ, Feng C. Outcome of small intestinal submucosa graft for repair of anterior urethral strictures. *Int J Urol*. 2013;20(6):622-629.

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

REFERANSLAR

1. Mantovani F, Tondelli E, Cozzi G, et al. Reconstructive urethroplasty using porcine acellular matrix (SIS): evolution of the grafting technique and results of 10-year experience. *Urologia*. 2011;78(2):92-97.
2. Palminteri E, Berdondini E, Fusco F, De Nunzio C, Salonia A. Long-term results of small intestinal submucosa graft in bulbar urethral reconstruction. *Urology*. 2012;79(3):695-701.
3. Xu YM, Fu Q, Sa YL, Zhang J, Song LJ, Feng C. Outcome of small intestinal submucosa graft for repair of anterior urethral strictures. *Int J Urol*. 2013;20(6): 622-629.

Keep dry

Manufacturer

Use-by date

Attention, see instructions for use

Do not re-use

Temperature limit

STERILE

EO

Sterilized using ethylene oxide

Magnetic resonance safe

Date of manufacture

MANUFACTURER

COOK BIOTECH
INCORPORATED
1425 Innovation Place
West Lafayette, IN 47906 U.S.A.

EC REPRESENTATIVE

Cook Biotech Europe ApS
Sandet 6, DK-4632
Bjaeverskov, DENMARK
Tel: +45 56868686

www.cookbiotech.eu
© COOK BIOTECH 2020

September 2020