**COOK®**  
MEDICAL

**CE** 0123

**EN**  
**3**

**OASIS® Extracellular Matrix**

Instructions for Use

**CS**  
**5**

**Extracelulární matrice OASIS®**

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**DA**  
**9**

**OASIS® ekstracellulær matrix**

Brugsanvisning

**DE**  
**12**

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**EL**  
**15**

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
## OASIS® EXTRACELLULAR MATRIX

### DEVICE DESCRIPTION

OASIS® EXTRACELLULAR MATRIX (ECM) is a dried sheet of small intestinal submucosa (SIS) used to manage partial-thickness wounds (comprising only epidermal and/or dermal layers) and full-thickness wounds (reaching subcutaneous tissue). Once applied, the natural composition of SIS allows for the patient's healing mechanisms to deposit cells and collagen during cellular and extracellular matrix metabolism. Typically within one week, the patient's tissues incorporate the SIS, resulting in new, remodeled tissue.<sup>1-3</sup> Results from a randomized controlled study showed that the mean frequency of OASIS reapplication was  $6.4 \pm 1.4$  days and demonstrated clinical efficacy with no adverse events at that frequency.<sup>3</sup>

### INTENDED USE

OASIS ECM is intended for the management of partial and full-thickness wounds. It is indicated to treat the following types of chronic and acute wounds:

- Abrasions
- Burns (second-degree only)
- Chronic vascular ulcers
- Diabetic ulcers
- Donor sites
- Lacerations
- Pressure ulcers
- Post-surgical dermal wounds
- Skin tears
- Venous ulcers

**Rx ONLY** This symbol means the following:

**CAUTION: This product is intended for use by trained medical professionals.**

**MR** This symbol means the following: Magnetic Resonance Safe

OASIS ECM is supplied sterile in peel-open packages and is intended for one-time use.

### COMPOSITION

OASIS ECM is derived using a process that retains the natural composition of matrix molecules such as collagen (Types I, III, IV, VI), glycosaminoglycans (hyaluronic acid, chondroitin sulfate A and B, heparin, and heparan sulfate), glycoproteins (fibronectin) and growth factors (FGF-2, TGF- $\beta$ ).<sup>4,5</sup>

### CONTRAINDICATIONS

OASIS ECM is derived from a porcine source and should not be used in patients with known sensitivity to porcine material. OASIS ECM is not indicated for use in third degree burns.

### PRECAUTIONS

- Do not re-sterilize. Discard all open and unused portions of the OASIS ECM.

- OASIS ECM is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- OASIS ECM must be used prior to the expiration date.
- Discard OASIS ECM if mishandling has caused possible damage or contamination.
- OASIS ECM should not be applied until excessive exudate, bleeding, acute swelling, and infection is controlled.
- No studies have been conducted to evaluate the performance of OASIS ECM in pediatric populations or in patients who are pregnant.

## POTENTIAL COMPLICATIONS

The following complications are possible. If any of these conditions occur, the device should be removed.

- Allergic reaction
- Blistering
- Infection
- Inflammation (Initial application of wound dressings may be associated with transient, mild, localized inflammation.)
- Pain
- Redness
- Swelling

## STORAGE

Store OASIS ECM in a clean, dry location at room temperature.

## STERILIZATION

OASIS ECM has been sterilized with ethylene oxide.

## INSTRUCTIONS FOR USE

**NOTE:** Always handle the OASIS ECM using aseptic technique.

1. Prepare the wound area using standard of care methods to ensure wound is free of debris and necrotic tissue. If necessary, surgically debride and treat the wound bed with standard of care techniques to ensure the wound shows granulation tissue and the edges contain viable tissue.
2. Trim the OASIS ECM to the appropriate size based on the size of wound area, ensuring that the graft covers the entire wound surface area and extends slightly beyond the wound margins. If the wound is larger than a single sheet, then multiple sheets may be used. Overlap adjoining sheets to provide coverage of the entire wound.
3. Apply the OASIS ECM in its dry state over the wound. Hydrate the device with sterile saline or lactated Ringer's solution. Alternatively, hydrate the OASIS ECM by placing it in a basin with sterile saline or lactated Ringer's solution prior to placing it over the wound. A hydration time of greater than 1 minute is not required.
4. Verify that the OASIS ECM extends beyond the margins of the wound to obtain the appropriate overlap. Smooth the OASIS ECM into place to ensure the sheet is in contact with the underlying wound bed. User may secure the graft with preferred method of fixation.
5. Apply a porous, non-adherent dressing over the OASIS ECM and affix using method of choice.

**NOTE:** A non-adherent dressing will help secure and protect the OASIS ECM during secondary dressing changes, while also

allowing fluid management to be controlled by secondary dressing layers.

6. Assess the wound to determine moisture conditions and apply appropriate secondary dressing layers to maintain a moist wound environment. Change the secondary dressing layers as needed to maintain a moist, clean wound area.

**NOTE:** To prevent damage to the newly-incorporated OASIS ECM, ensure that the non-adherent dressing (step 5) is only changed when ready to assess the wound.

7. Discard any unused portions of unused OASIS ECM according to institutional guidelines for biological waste. Reapply, as needed, if the OASIS ECM is no longer visible. Typically, reapplication is needed every 3-7 days until the wound is re-epithelialized.
8. As healing occurs, sections of the OASIS ECM may gradually peel. Carefully remove any remaining loose product around the edge as needed.

**IMPORTANT:** If a gel forms on the wound surface, do not attempt to forcibly remove it. Successful integration of the OASIS ECM may form a caramel-colored or off-white gel. Do not remove this gel by debridement. This caramelization contains ECM, which continues to replace deficient and missing ECM in the wound.

9. Gently cleanse the wound surface with sterile saline; leave the ECM gel intact.
10. Carefully reassess the wound and record healing progression and other relevant information.

## REFERENCES

1. Shi L, Ronfard V. Biochemical and biomechanical characterization of porcine small intestinal submucosa (SIS): a mini review. *Int J Burn Trauma*. 2013;3:173-179.
2. Yeh DD, Nazarian RM, Demetri L, et al. Histopathological assessment of OASIS Ultra on critical-sized wound healing: a pilot study. *J Cutan Pathol*. 2017;44(6):523-529.
3. Romanelli M, Dini V, Bertone M, Barbanera S, Brilli C. OASIS wound matrix versus Hyaloskin in the treatment of difficult-to-heal wounds of mixed arterial/venous etiology. *Int Wound J*. 2007;4(1):3-7.
4. Hodde J, Janis A, Ernst D, Zopf D, Sherman D, Johnson C. Effects of sterilization on an extracellular matrix scaffold: Part I. Composition and matrix architecture. *J Mater Sci Mater Med*. 2007;18(4):537-543.
5. Hodde JP, Badylak SF, Brightman AO, Voytik-Harbin SL. Glycosaminoglycan content of small intestinal submucosa: A bioscaffold for tissue replacement. *Tissue Eng*. 1996;2(3): 209-217.

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## LITERATURA

1. Shi L, Ronford V. Biochemical and biomechanical characterization of porcine small intestinal submucosa (SIS): a mini review. *Int J Burn Trauma.* 2013;3:173-179.
2. Yeh DD, Nazarian RM, Demetri L, et al. Histopathological assessment of OASIS Ultra on critical-sized wound healing: a pilot study. *J Cutan Pathol.* 2017;44(6):523-529.
3. Romanelli M, Dini V, Bertone M, Barbanera S, Brilli C. OASIS wound matrix versus Hyaloskin in the treatment of difficult-to-heal wounds of mixed arterial/venous etiology. *Int Wound J.* 2007;4(1):3-7.
4. Hodde J, Janis A, Ernst D, Zopf D, Sherman D, Johnson C. Effects of sterilization on an extracellular matrix scaffold: Part I. Composition and matrix architecture. *J Mater Sci Mater Med.* 2007;18(4):537-543.
5. Hodde JP, Badylak SF, Brightman AO, Voytik-Harbin SL. Glycosaminoglycan content of small intestinal submucosa: A bioscaffold for tissue replacement. *Tissue Eng.* 1996;2(3): 209-217.

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2. Yeh DD, Nazarian RM, Demetri L, et al. Histopathological assessment of OASIS Ultra on critical-sized wound healing: a pilot study. *J Cutan Pathol.* 2017;44(6):523-529.
3. Romanelli M, Dini V, Bertone M, Barbanera S, Brilli C. OASIS wound matrix versus Hyaloskin in the treatment of difficult-to-heal wounds of mixed arterial/venous etiology. *Int Wound J.* 2007;4(1):3-7.
4. Hodde J, Janis A, Ernst D, Zopf D, Sherman D, Johnson C. Effects of sterilization on an extracellular matrix scaffold: Part I. Composition and matrix architecture. *J Mater Sci Mater Med.* 2007;18(4):537-543.
5. Hodde JP, Badylak SF, Brightman AO, Voytik-Harbin SL. Glycosaminoglycan content of small intestinal submucosa: A bioscaffold for tissue replacement. *Tissue Eng.* 1996;2(3): 209-217.

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## REFERÊNCIAS

1. Shi L, Ronford V. Biochemical and biomechanical characterization of porcine small intestinal submucosa (SIS): a mini review. *Int J Burn Trauma*. 2013;3:173-179.
2. Yeh DD, Nazarian RM, Demetri L, et al. Histopathological assessment of OASIS Ultra on critical-sized wound healing: a pilot study. *J Cutan Pathol*. 2017;44(6):523-529.
3. Romanelli M, Dini V, Bertone M, Barbanera S, Brilli C. OASIS wound matrix versus Hyaloskin in the treatment of difficult-to-heal wounds of mixed arterial/venous etiology. *Int Wound J*. 2007;4(1):3-7.
4. Hodde J, Janis A, Ernst D, Zopf D, Sherman D, Johnson C. Effects of sterilization on an extracellular matrix scaffold: Part I. Composition and matrix architecture. *J Mater Sci Mater Med*. 2007;18(4):537-543.

5. Hodde JP, Badylak SF, Brightman AO, Voytik-Harbin SL.  
Glycosaminoglycan content of small intestinal submucosa:  
A bioscaffold for tissue replacement. *Tissue Eng.* 1996;2(3):  
209-217.

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## REFERENSER

1. Shi L, Ronford V. Biochemical and biomechanical characterization of porcine small intestinal submucosa (SIS): a mini review. *Int J Burn Trauma*. 2013;3:173-179.
2. Yeh DD, Nazarian RM, Demetri L, et al. Histopathological assessment of OASIS Ultra on critical-sized wound healing: a pilot study. *J Cutan Pathol*. 2017;44(6):523-529.
3. Romanelli M, Dini V, Bertone M, Barbanera S, Brilli C. OASIS wound matrix versus Hyaloskin in the treatment of difficult-to-heal wounds of mixed arterial/venous etiology. *Int Wound J*. 2007;4(1):3-7.
4. Hodde J, Janis A, Ernst D, Zopf D, Sherman D, Johnson C. Effects of sterilization on an extracellular matrix scaffold: Part I. Composition and matrix architecture. *J Mater Sci Mater Med*. 2007;18(4):537-543.
5. Hodde JP, Badylak SF, Brightman AO, Voytik-Harbin SL. Glycosaminoglycan content of small intestinal submucosa: A bioscaffold for tissue replacement. *Tissue Eng*. 1996;2(3): 209-217.

**Manufacturer**

**Use-by date**

**Attention, see instructions for use**

**Do not re-use**

**Temperature limit**

**Magnetic resonance safe**

**STERILE**

**EO**

**Sterilized using  
ethylene oxide**

**Date of manufacture**

**MANUFACTURER**

COOK BIOTECH  
INCORPORATED  
1425 Innovation Place  
West Lafayette, IN 47906 U.S.A.

**EC REPRESENTATIVE**

Cook Biotech Europe ApS  
Sandet 6, DK-4632  
Bjaeverskov, DENMARK  
Tel: +45 56868686