
COOK
MEDICAL

CE 0123

EN 3	Biodesign® Hernia Graft Instructions for Use
CS 5	Kýlní implantát Biodesign® Návod k použití
DA 7	Biodesign® hernieimplantat Brugsanvisning
DE 9	Biodesign® Hernienimplantat Gebrauchsanweisung
EL 11	Μόσχευμα κήλης Biodesign® Οδηγίες χρήσης
ES 13	Biodesign® Injerto para hernias Instrucciones de uso
FR 15	Implant biologique de hernie Biodesign® Mode d'emploi
HU 18	Biodesign® sérvgraft Használati utasítás
IT 20	Innesto per ernia Biodesign® Istruzioni per l'uso
NL 22	Biodesign® herniaprothese Gebruiksaanwijzing
NO 24	Biodesign® brokkimplantat Bruksanvisning
PL 26	Wszczep do przepuklin Biodesign® Instrukcja użycia
PT-BR 28	Enxerto de Hérnia Biodesign® Instruções de uso
SV 31	Biodesign® bräckgraft Bruksanvisning

FP0117-01B

BIODESIGN® HERNIA GRAFT

DEVICE DESCRIPTION

Biodesign® Hernia Graft is a multi-layered sheet of extracellular matrix derived from the small intestinal submucosa (SIS) of pigs. SIS is obtained from intestine using a process that retains the natural composition of matrix molecules.^{1,2}

Once implanted, the natural composition of SIS allows for the patient's healing mechanisms to deposit cells and collagen during cellular and extracellular matrix metabolism. Over the course of several months, the patient's tissues completely incorporate the SIS until it can no longer be detected, resulting in new remodeled tissue.³

INTENDED USE

The Cook® Biodesign® Hernia Graft is intended for implantation to reinforce soft tissues where weakness exists during ventral hernia repair. The graft is supplied sterile and is intended for one-time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

MR This symbol means the following: Magnetic Resonance Safe

HERNIA GRAFT This symbol means the following: Hernia Graft

This graft is intended for use by trained medical professionals.

CONTRAINDICATIONS

This graft is derived from a porcine source and should not be used in patients with known sensitivity to porcine material.

PRECAUTIONS

- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- **Do not resterilize.** Discard all open and unused portions of the graft.
- The graft is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard graft if mishandling has caused possible damage or contamination, or if the graft is past its expiration date.
- Ensure that graft is rehydrated prior to cutting, suturing, stapling, tacking or loading of the graft laparoscopically.
- Ensure that all layers of the graft are secured when suturing, stapling, or tacking.
- Place graft in maximum possible contact with healthy, well-vascularized tissue to encourage cell ingrowth and tissue remodeling.
- Suturing, stapling, or tacking more than one graft together may decrease graft performance.
- No studies have been conducted to evaluate the reproductive impact of the clinical use of the graft.
- Extended rehydration or excessive handling could lead to partial delamination of superficial layers of the graft.
- Care should be taken when device is placed in infected wounds.
- Care should be taken to avoid damage to the graft when loading laparoscopically. It is recommended to load through a 10 mm or larger port.
- If wound is left open, keep graft moist to prevent dryness.

POTENTIAL COMPLICATIONS

Possible adverse reactions with the use of biologic hernia prostheses may include:

- adhesion
- allergic reaction
- bowel erosion
- bowel obstruction
- discharge
- fever
- fistula formation
- hematoma
- hernia recurrence
- infection/abscess formation
- inflammation
- pain
- premature degradation
- seroma formation
- wound dehiscence

NOTE: Peer-reviewed published literature suggests that patients with innate collagen deficiencies and recurrent hernias are more susceptible to subsequent hernia formation.⁴

NOTE: Peer-reviewed published literature demonstrates that the risk of complications, such as delayed wound infection, premature degradation, hernia recurrence, bowel erosion, and the need for re-operation, are higher in patients who are critically ill or who have severely contaminated abdomens when they are treated with biologic hernia grafts in the presence of infection.⁵ In patients who have severely contaminated abdomens, a staged approach to

hernia repair in which the infection is controlled before definitive abdominal closure is attempted may result in improved outcomes.

STORAGE

This graft should be stored in a clean, dry location at room temperature.

STERILIZATION

This graft has been sterilized with ethylene oxide.

INSTRUCTIONS FOR USE

Required Materials

- Sterile forceps
- Rehydration fluid: room temperature, sterile saline or sterile lactated Ringer's solution.

NOTE: Always handle the graft using aseptic technique, minimizing contact with latex gloves.

1. Using aseptic technique, remove the sterile inner package from the outer packaging and place it in the sterile field.
2. Open the sterile inner package carefully. Aseptically, remove the graft from the inner package.
3. Rehydrate the graft in room temperature, sterile saline or sterile lactated Ringer's solution until the desired handling characteristics are achieved. For the Hernia Graft, rehydrate for less than one (1) minute prior to cutting, suturing, stapling, tacking or loading the device laparoscopically. If the graft is provided in a tray, it can be rehydrated in the tray.
4. Prepare the hernia repair site using standard surgical techniques.

MESH PLACEMENT

- Preperitoneal/retrorectus (Rives and Stoppa technique) or intraperitoneal placement (underlay) is recommended for open procedures.
 - Intraperitoneal placement (underlay) is recommended for a laparoscopic approach.
5. Using aseptic technique, trim the graft to fit the site, providing an allowance for overlap. Fundamental surgical principles suggest recurrence can be minimized if the mesh overlaps the surrounding tissue by at least 4-5 cm in all directions.⁶

NOTE: If the graft is cut too small for the defect, excess tension may be placed on the suture line. This can result in recurrence of the original hernia or development of a defect in the adjacent tissues.

6. Using aseptic technique, transfer the graft to the surgical site and suture, staple, or tack into place, avoiding excess tension.

NOTE: Excessive handling may lead to partial delamination of a superficial layer. This can be managed by securing the detached portion to adjacent tissue with sutures or staples.

NOTE: Surgical experience indicates that suturing or stapling the graft with close tissue approximation produces better outcomes. Leave less than or equal to 3cm between sutures with bite depth of 1cm. Use permanent or long-term absorbable sutures or tacks according to surgeon preference.

NOTE: Interrupted sutures can provide additional security against recurrence of tissue defect in the event of suture failure.

NOTE: The Hernia Graft has been designed with an open matrix structure to optimize tissue incorporation and angiogenesis.

NOTE: It is recommended to minimize the bridging of defects to improve outcomes.

7. Complete the standard surgical procedure.
8. Discard any unused portions according to institutional guidelines for medical waste.
9. Place closed suction drains for 2-6 weeks. Remove when output is less than 20 mL/24 hours for at least two (2) consecutive days or until drain is dry.

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West Lafayette, IN 47906 U.S.A.
AUSTRALIA PTY. LTD.