
COOK
MEDICAL

CE 0123

- EN**
3 **Biodesign®**
Hiatal Hernia Graft
Instructions for Use
- CS**
4 **Biodesign®**
Graft pro hiátového hernii
Návod k použití
- DA**
6 **Biodesign®**
Hiatushernieimplantat
Brugsanvisning
- DE**
8 **Biodesign®**
Hiatushernien-Implantat
Gebrauchsanweisung
- EL**
10 **Biodesign®**
Μόσχευμα διαφραγματοκήλης
Οδηγίες χρήσης
- ES**
12 **Biodesign®**
Injerto para hernias de hiato
Instrucciones de uso
- FI**
13 **Biodesign®**
Hiatusherniasiiirre
Käyttöohjeet
- FR**
15 **Biodesign®**
Implant biologique pour hernie hiatale
Mode d'emploi
- IT**
17 **Biodesign®**
Innesto per ernia iatale
Istruzioni per l'uso
- NL**
19 **Biodesign®**
Hiatusherniaprothese
Gebruiksaanwijzing
- NO**
21 **Biodesign®**
Spiserørsbrokkimplantat
Bruksanvisning
- PT-BR**
23 **Biodesign®**
Enxerto para hérnia de hiato
Instruções de uso
- SV**
25 **Biodesign®**
Hiatusbräckgraft
Bruksanvisning
- TR**
26 **Biodesign®**
Hiatal Herni Grefti
Kullanma Talimatı

FP0141-01B

BIODESIGN® HIATAL HERNIA GRAFT

DESCRIPTION

The Biodesign® Hiatal Hernia Graft is a multi-layered sheet of extracellular collagen matrix derived from porcine small intestinal submucosa (SIS). The Biodesign Hiatal Hernia Graft is implanted by a trained medical professional to support weakened soft tissue during open or laparoscopic surgery to repair hiatal hernias. Once implanted, the natural composition of SIS allows for the patient's healing mechanisms to deposit cells and collagen during cellular and extracellular matrix metabolism. Over the course of 7 to 17 months, depending on the implant location, the patient's tissues completely incorporate the SIS resulting in new, remodeled tissue.¹⁻³

INTENDED USE

The Biodesign Hiatal Hernia Graft is used for implantation to reinforce soft tissues where weakness exists, during hiatal hernia repair. This graft is supplied sterile in peel-open packages and is intended for one-time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

MR This symbol means the following: Magnetic Resonance Safe

HIATAL HERNIA GRAFT This symbol means the following: Hiatal Hernia Graft

This product is intended for use by trained medical professionals.

CONTRAINDICATIONS

The graft is derived from a porcine source and should not be used for patients with known sensitivity to porcine material.

PRECAUTIONS

- The graft is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to graft failure and/or transmission of disease.
- **Do not resterilize.** Discard all open and unused portions.
- The graft is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard the graft if mishandling has caused possible damage or contamination, or if the graft is past its expiration date.
- Ensure that the graft is hydrated prior to cutting or suturing. Insufficient hydration may cause graft damage during laparoscopic delivery.
- Ensure that all layers of the graft are secured when suturing.
- Ensure that the graft does not touch the esophagus.

POTENTIAL COMPLICATIONS

The following complications are possible with the use of surgical graft materials:

- Adhesion
- Allergic reaction
- Fistula formation
- Graft detachment/migration
- Hematoma
- Infection
- Inflammation
- Recurrence of tissue defect
- Seroma formation

The following complications are associated with use of graft materials in hiatal hernia repair:

- Dysphagia
- Esophageal erosion
- Perforation
- Stricture formation

STORAGE

Store the graft in a clean, dry location at room temperature.

STERILIZATION

The graft has been sterilized with ethylene oxide.

INSTRUCTIONS FOR USE

Required Materials

- Sterile basin
- Sterile smooth forceps
- Hydration fluid: room temperature sterile saline or sterile lactated Ringer's solution

NOTE: Handle the graft using aseptic technique, minimizing contact with latex gloves.

1. Remove the packaging containing the graft from the envelope.
2. Using aseptic technique, open the outer pouch and pass the inner pouch containing the graft onto the sterile field.
3. Place the graft into the sterile basin on the sterile field. Multiple grafts may be hydrated simultaneously in the same basin.
4. Place enough hydration fluid in the basin to completely cover the graft(s).
5. Hydrate the graft in the fluid until the desired handling characteristics are achieved. A hydration time of greater than one minute is not required.
6. Prepare the graft site using standard surgical techniques.

NOTE: Suturing the graft with close tissue approximation produces better outcomes. Maximum graft-tissue approximation is only achieved with complete crural closure.

7. If necessary, trim the graft to fit the patient's anatomy, providing a small allowance for overlap with the adjacent tissues.


NOTE: The graft may be trimmed prior to hydration to fit the patient's anatomy. Be sure to hydrate the graft prior to suturing it into place. See step 5.

NOTE: For laparoscopic placement, hydrate the graft prior to loading it into the delivery port. See step 5.

8. Suture the graft into place, avoiding excess tension.

NOTE: Clinical literature suggests the use of tacks or staples at the diaphragm may be associated with cardiovascular complications.⁴

NOTE: The 7X10-U graft has a larger hole in the corner to identify the graft orientation. The larger hole should be placed at the patient's upper left.

9. Complete the surgical procedure.
10. Discard unused portions of the graft according to institutional guidelines for disposal of medical waste.

REFERENCES

1. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol.* 2004;172(1):215-218.
2. Franklin ME, Jr., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinal submucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia.* 2004;8(3):186-189.
3. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestine submucosa as a graft material for pubovaginal slings: durability and results. *Urology.* 2003;62(5):805-809.
4. Köckerling F, Schug-Pass C, Bittner R. A word of caution: never use tacks for mesh fixation to the diaphragm! *Surg Endosc.* 2018;32(7):3295-3302.

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LITTERATUR

1. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol.* 2004;172(1):215-218.
2. Franklin ME, Jr., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinal submucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia.* 2004;8(3):186-189.
3. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestine submucosa as a graft material for pubovaginal slings: durability and results. *Urology.* 2003;62(5):805-809.
4. Kockerling F, Schug-Pass C, Bittner R. A word of caution: never use tacks for mesh fixation to the diaphragm! *Surg Endosc.* 2018;32(7):3295-3302.

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QUELLEN

1. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol.* 2004;172(1):215-218.
2. Franklin ME, Jr., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinal submucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia.* 2004;8(3):186-189.

3. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestine submucosa as a graft material for pubovaginal slings: durability and results. *Urology*. 2003;62(5):805-809.

4. Köckerling F, Schug-Pass C, Bittner R. A word of caution: never use tacks for mesh fixation to the diaphragm! *Surg Endosc*. 2018;32(7):3295-3302.

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REFERANSER

1. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol.* 2004;172(1):215-218.

2. Franklin ME, Jr., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinal submucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia*. 2004;8(3):186-189.
3. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestine submucosa as a graft material for pubovaginal slings: durability and results. *Urology*. 2003;62(5):805-809.
4. Köckerling F, Schug-Pass C, Bittner R. A word of caution: never use tacks for mesh fixation to the diaphragm! *Surg Endosc*. 2018;32(7):3295-3302.

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


REFERANSLAR

1. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol.* 2004;172(1):215-218.
2. Franklin ME, Jr, Gonzalez JJ, Jr, Glass JL. Use of porcine small intestinal submucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia.* 2004;8(3):186-189.
3. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestine submucosa as a graft material for pubovaginal slings: durability and results. *Urology.* 2003;62(5):805-809.
4. Köckerling F, Schug-Pass C, Bittner R. A word of caution: never use tacks for mesh fixation to the diaphragm! *Surg Endosc.* 2018;32(7):3295-3302.

Magnetic resonance safe

Keep dry

Manufacturer

Use-by date

Attention, see instructions for use

Do not re-use

Temperature limit

Sterilized using ethylene oxide

Date of manufacture

MANUFACTURER

COOK BIOTECH
INCORPORATED
1425 Innovation Place
West Lafayette, IN 47906 U.S.A.

EC REPRESENTATIVE

Cook Biotech Europe ApS
Sandet 6, DK-4632
Bjaeverskov, DENMARK
Tel: +45 56868686